

150A Biochemical Laboratory Incubator, Biochemical Chamber



Basic Info

- Model NO.: 150A
- Structure: Floor Type
- Certification: ISO, SGS
- Type: Incubator
- Specification: 1pc in 1 box

- Customized: Customized
- Material: Steel
- Application: School, Hospital, Lab
- Trademark: JY/OEM
- Origin: Changzhou, China

Product Description

Biochemical Laboratory Incubator

1. Easy to use, steady performance, good durability
2. ISO13485/14001/9001 Certified
3. Since 1984

Biochemical incubator is widely used for storage and culture of bacteria, mold, microorganism, tissue and cell, water quality analysis and BOD test. Applicable for biology, genetic engineering, medicine, hygiene and disease control, environmental protection and other sectors.

Datasheet:

Specifications of 150A Intelligent Biochemical Laboratory Incubator	
Power Supply:	AC 220V/110V 50Hz/60Hz
Volume:	150L
Temperature Range:	5-50celsius
Temperature fluctuation:	± 0.5celsius
Lighting function:	Yes
Digital Display:	Yes
Widely applicable for lab, biology, chemistry, physical, hospital, food and environmental protection fields.	