

501 Ultra Thermostatic Water Bath



Basic Info

- Model NO.: 501
- Structure: Desktop
- Certification: CE, ISO, SGS
- Type: Heating Type
- Trademark: JY
- Origin: Changzhou, China
- Customized: Customized
- Material: Stainless Steel
- Application: School, Hospital, Lab
- 350X240X180mm: Outer Size:420X350X280mm
- Specification: 1PC in 1 carton
- HS Code: 8419899090

Product Description

Laboratory Water Bath

1. S steel inner, LED, Thermostatic
2. ISO13485/14001/9001 Certified
3. Since 1984

We are focused on design and manufacturing laboratory device and related products over 20 years. We export to countries all over the world like the middle east and south Africa.... with the most competitive prices and quality products. Except for standard, we are also welcome customization according to your design and style.

Specifications of 501 Ultra Thermostatic Lab Water Bath	
Power Supply	AC 220V/110V 50Hz/60Hz
Temperature control range	Room temperature to 95 Celsius
Water temp uniformity	≤ 0.5celsius
Working room	35X24X18cm
Pump capacity	6L/min
Digital display	Yes
Easy to use, rich function configuration, high quality, wide applicability. Stainless steel made water tank and pump. OEM service is offered.	