

601 Ultra Thermostatic Digital Laboratory Water Bath



Basic Info

- Model NO.: 601
- Structure: Desktop
- Application: School, Hospital, Lab
- Super Water Bath: 8 Liter
- Specification: CE, ISO13485/14001/9001 Certified Factory
- HS Code: 8419899090
- Customized: Customized
- Certification: CE, ISO, SGS
- Type: Heating Type
- Trademark: JY/OEM
- Origin: China

Product Description

601 Ultra thermostatic Digital Laboratory Water Bath

1. SUS inner and lid, pump capacity: 6L/h, precise, digital display
2. CE, ISO13485/14001/9001 Certified Factory

Specifications of 601 Ultra thermostatic Digital Laboratory Water Bath	
Power Supply:	AC 220V/110V 50Hz/60Hz
Temperature control range:	Room temperature+5 to 100 Celsius
Fluctuate of Constant Temperature:	± 0.1 Celsius
Digital temperature control:	Yes
Power:	1000W
Pump Flow:	6L/min
Working room size:	L28 W24 H18cm
Easy to use, rich function configuration, high quality, wide applicability. Stainless steel made water tank. OEM service is offered.	