

76-1A Digital Display Glass Type Thermostatic Water Bath



Basic Info

- Model NO.: 76-1A
- Structure: Desktop
- Certification: ISO, SGS
- Type: Heating Type
- Specification: 1PC in 1 carton

- Customized: Customized
- Material: Glass
- Application: School, Hospital, Lab
- Trademark: JY
- Origin: Changzhou, China

Product Description

Glass Type Thermostatic Water Bath

1. RT-100Celsius, LED temp control, Glass type
2. ISO13485/14001/9001 Certified
3. Since 1984

Specifications of 76-1A Glass Type Thermostatic Water Bath

Power Supply	AC 220V/110V 50Hz/60Hz
Temperature control range	Room temperature to 100 Celsius
Heating power	1000W
Motor power	60W
Speed	0-2000rpm
Digital display	Temp control
Easy to use, rich function configuration, high quality, wide applicability. Glass tank. OEM service is offered.	