
79-2 Laboratory Bi-Directional Magnetic Heated Stirrer



Basic Info

Model NO.: 79-2Trademark: JY

• Origin: Changzhou, China

Type: Magnetic StirrerSpecification: 1 pc in 1 boxHS Code: 8479820090

Product Description

78-1 Laboratory Magnetic Stirrer With Hotplate/ Bio-mixer

1. Simple & easy to use, reliable, bi-directional

2. ISO13485/14001/9001

3. Since 1984

Magnetic Hotplate Stirrer, an ideal essential device for laboratories, is widely applied in colleges and universities, environment protecting, researching, sanitation, antilipidemic, petroleum, metallurgy, chemical and medical units.

Specifications of 79-2 Laboratory Bi-directional Magnetic Heated Stirrer	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	0-2400rpm
Motor powers	40W
Heating power	300W