## Centrifuge

## Outline of uses:

80-2Bisan improvementon80-1A.It is possessed of these special feature following: it is beautiful in appearance, large capacity but small in bulk and has a complete function and a stable performance, You can chose the speed that you need at the set range, And it can keep balance automaticly, in addition, it has a little rate of increase in temperature and can be used efficiently and widely, It is a dependable

80-2B Centrifuge

instrument for qualitative analysis of serum, plasma and immunity in hospital sorchemical and biochemical laboratories.

## Specification:

- 1. PowerSupply:AC220V 50Hz
- 2. Speed:0Ђ4000rpm
- 3. Timer range:0Ђ60min
- 4. Max.RCF:2325xg
- 5. Capacity:12¦20ml
- 6. Motorpower:40 w
- 7. Weight:9.5Kg