

800 Series Electric Centrifuge

Centrifuge

Features:

Electric centrifugal machine has features of artistic shape, large capacity, small volume and complete function. It has such advantages as stable performance, adjustable speed and automatic adjustment and balance, lower temperature, high use efficiency and extensive applicability.

This product applies into the qualitative analysis of blood serum, blood plasma and car amide for medical product, blood station, and clinic test as well as biochemical lab.



800Centrifuge



800DTypeElectricCentrifugal



80-1TypeElectricCentrifugal



80-2TypeCentrifugalSedimentator

Specification:

Type	Description	RotationSpeed	Capacity	MaximumRelativeDisplacement	PowerSupply
800&800D	ElectricCentrifugalMachine	0-4000rpm	15mlx6	1430xg	220V50HZ (110V60Hzoption)
80-1	ElectricCentrifugalMachine	0-4000rpm	20mlx6	1430xg	
80-2	CentrifugalSedimentator	0-4000rpm	20mlx12	2325xg	

Remark:800Dwithtimefunction