

## 85-2 Laboratory Magnetic Hotplate Stirrer



### Basic Info

- Model NO.: 85-2
- Stainless Steel Plate: with Stirrer
- Specification: 1 pc in 1 box
- HS Code: 8479820090
  
- Type: Magnetic Stirrer
- Trademark: JY
- Origin: China

### Product Description

85-2 Laboratory Magnetic Stirrer with Hotplate

1. Simple&easy use, reliable, large capacity, LED
2. ISO13485/14001/9001
3. Since 1984

Magnetic Hotplate Stirrer, an ideal essential device for laboratories, is widely applied in colleges and universities, environment protecting, researching, sanitation, antilipidemic, petroleum, metallurgy, chemical and medical units.

Specifications of 85-2 Digital Laboratory Magnetic Stirrer	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	0-2400rpm
Motor powers	40W
Heating power	250W
Temperature control range	Room temperature to 150 Celsius
Timer	0-120min
Two-way rotation	Yes
LED display	Yes, temperature and speed indication

capacity	20-3000ml
Easy to use, rich function configuration, cost-efficient.	