

85-2b Heating Magnetic Stirrer for Laboratory



Basic Info

- Model NO.: 85-2B
- Type: Magnetic Stirrer
- Trademark: JY
- Origin: China
- Application: Lab
- Easy to Operate: Digital Display
- Specification: ISO13485/14001/9001
- HS Code: 8479820090

Product Description

85-2B Laboratory Magnetic Stirrer With Hotplate/Biomixer

1. LED, wide applicability, large capacity, timer, speed and temp digital display
2. ISO13485/14001/9001
3. Since 1984

Magnetic Hotplate Stirrer, an ideal essential device for laboratories, is widely applied in colleges and universities, environment protecting, researching, sanitation, antilipidemic, petroleum, metallurgy, chemical and medical units.

Specifications of 85-2B Laboratory Magnetic Hotplate Stirrer	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	0-2400rpm
Motor powers	40W
Heating power	300W
Temperature control range	Room temperature to 100 celsius
Timer	0-120min
Two-way rotation	Yes

LED display	Yes, temperature and speed indication
Easy to use, rich function configuration, cost-efficient.	