90-2 Centrifuge(Digital)

Outline of uses:

90-2 is an improvement on 80-2B. It is possessed of these special feature following: It is beautiful in appearance, large capacity but small in bulk and has a complete function and a stable performance, you can chose the speed that you need at the set range, and it can keep balance automatically, In addition, it has a little rate of increase in temperature and can be used efficiently and widely, It is a

dependable instrument for qualitative analysis of serum, plasma and immunity in hospitals or chemical and biochemical laboratories.

	Model	MaxSpeed	Capacity	MaxRCF	TimerRange	NetWeight
	90-2	4000r/min	20ml*12	2325g	0~60min	9Kg
ı						