

ASE-24 (Round)



Features:

Each independent valve control durable, more convenient operation

Anti-cross-contamination, anti-fog glass vacuum tank

Can be equipped with the large capacity acquisition, and rapid enrichment of the drying unit, batch processing of samples.

We recommend that you buy solid-phase extraction, matching the AP series I produced oil-free vacuum pump.

ASE-24

Pressure capability: better than 80kpa

Can handle the number of samples: 24

Supporting test tube rack size: $\Phi 11$, $\Phi 13$, $\Phi 15$

Length of 105mm and below. ">