Electrical Heating Incubator

----Intelligentized Microprocessor Control

---LCD Display

Features:

- 1. With mirror surface in stainless steel, four square semicircle transition, shelf support handling freely, easy for cleaning within the case.
- 2. Giant-screen liquid crystal display, groups-data display in one screen, special all directions navigation key, easily operation, with more human nature in design.
- ${\it 3.}\ Use\ new\ handling\ of\ fuzzy\ logic\ controller,\ fast\ to\ set value,\ accuracy\ and\ stable\ in\ operation.$
- 4. Independent temperature-limited alarm system can automatically interrupt when it is over limit temperature, ensuring the safe operation of the experiment, with no accident. (Option)
- 5. RS485 interface may link computer and remote monitor the condition of machine.(Option)

Specification:

Model	BPH-9042	BPH-9082	BPH-9162	BPH-9272
Electrical Requirements	220V 50Hz			
Temperature Range	RT+5 - 80 °C			
Temperature Control	0.1 °C			
Temperature Stability	□,1 °C			
Ambient Temperature	+5 ~ 40 °C			
Input Power	250W	350W	550W	750W
Interior Dimension	350*320*350	400*400*500	500*400*750	520*460*1050
(mm)W*D*H				
Shelves per chamber(Normal)	2 PCS			
Timing Range	1~9999min			

