

BWS-12 Series



Precise Water Bath

Features

- * Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users of time.
- * LCD screen, multiple data display with timing function, easy to operate.
- * With timer function.
- * A stamping molding stainless steel tank, easy to clean.
- * Steepletop stainless steel cover to avoid water evaporates and water enters the sample.

Option:

- * Printer
- * RS485 connector

Specification:

Model	BWS-12	BWS-12G (With electromagnetic pump)	BWS-27	BWS-27G (With electromagnetic pump)
Electrical Requirements	220V 50HZ			
Power Consumption	800W		1000W	
Temperature Range	RT+5~80			
Temperature stability	0.2			
Tracking Alarm	0.1			
Chamber Volume	11L		20L	
Chamber Dimension (W*D*H, mm)	300*240*200		500*300*200	
Timing Range	1~999min			