

BWS-5 Series



Features

- * A stamping molding stainless steel tank, easy to clean.
- * LCD screen, multiple data display with timing function, easy to operate.
- * Stainless steel shelves cover heater and sensor to avoid damage during use.
- * Once-forming stainless steel lid.
- * Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time.
- * Independent temperature-limiting alarm system
- * Temperature error alarm
- * Test tube holder can be placed. (Optional)

Spec:

Model	BWS-5	BWS-10	BWS-20	BWS-0505	BWS-0510
Electrical Requirement	AC220V 50HZ				
Power Consumption	500W	1000W	1500W	500W+500W	500W+1000W
Temperature Range	RT+5~99°C				
Temperature Stability	±0.5°C				
Temp Alarm	±2°C				
Interior Dimension (W×D×H, mm)	130*280*150	220*280*150	290*490*150	130*280*150 130*280*150	130*280*150 290*490*150
Timing Range	1~9999min				
Remark	2holes Φ12mm	4holes Φ92mm	6holes Φ92mm	2holes +2holes	2holes +4holes

Option:

Intelligent Programmable temperature controller
 Printer
 Independent temperature-limiting Alarm system
 RS485 connector