

## **BXM-30R**



Product Name: VERTICAL AUTOCLAVE

Product Model: BXM-30R

### APPLICABLE SCOPE:

The Autoclave is the equipment which utilizes the pressure saturated steam for sterilizing goods quickly and reliable. It is applicable to be used in medical and health work, scientific and agriculture institutes. It is a perfect equipment for sterilizing the medical equipments, dressings, glassware, liquor culture mediums, etc.

### PRODUCT CHARACTERISTICS:

Self-control type

Outer case, inner vessel and baskets are made of SUS304, resistant to acid, alkali and corrosion

Microcomputer intelligent control

Over temperature protection device

Self expanding type seal ring, exhaust cold air automatically

Low water level alarm, self-control when lack of water

Self-exhaust when over pressure

Buzzer reminds when sterilization ends

### TECHNICAL PARAMETERS:

Model: BXM-30R

Volume: 30L

Power: 2KW

Power Supply: 220V±10% 50Hz±2%

Max. Working/Design Temp.: 126°C/128°C

Max. Working/Design Pressure: 0.14MPa/0.165MPa

Time Range (min): 0-120

Chamber Size (mm): Φ300×400

Basket Size (mm): Φ280×360×1pc

Overall Size (mm): 404×404×810