





Cell-02							
Picture and Name	Outside Dim H*W*D(mm)	Path Length	Inside Width	Volume (ml)	Product Number		
		(mm)	(mm)		ES Quartz Glass 190nm-2500nm	IR Quartz Glass 250nm-3500nm	Optical Glass 340nm-2500nm
 Micro cell with black walls and lid	45*12.5*7.5	5	2	0.35	Q-43	I-43	G-43
	45*12.5*12.5	10	2	0.7	Q-44	I-44	G-44
	45*12.5*22.5	20	2	1.4	Q-45	I-45	G-45
	45*12.5*32.5	30	2	2.1	Q-46	I-46	G-46
	45*12.5*42.5	40	2	2.8	Q-47	I-47	G-47
	45*12.5*52.5	50	2	3.5	Q-48	I-48	G-48
	45*12.5*102.5	100	2	7	Q-49	I-49	G-49
 Micro cell with black walls and teflon stopper	45*12.5*7.5	5	2	0.35	Q-53	I-53	G-53
	45*12.5*12.5	10	2	0.7	Q-54	I-54	G-54
	45*12.5*22.5	20	2	1.4	Q-55	I-55	G-55
	45*12.5*32.5	30	2	2.1	Q-56	I-56	G-56
	45*12.5*42.5	40	2	2.8	Q-57	I-57	G-57
	45*12.5*52.5	50	2	3.5	Q-58	I-58	G-58
	45*12.5*102.5	100	2	7	Q-59	I-59	G-59
 Self masking continuous flowthrough cell	45*12.5*7.5	5	φ3	0.035	Q-63	I-63	G-63
	45*12.5*12.5	10	φ3	0.07	Q-64	I-64	G-64
	45*12.5*22.5	20	φ3	0.14	Q-65	I-65	G-65
	45*12.5*32.5	30	φ3	0.21	Q-66	I-66	G-66
	45*12.5*7.5	5	φ2	0.015	Q-67	I-67	G-67
	45*12.5*12.5	10	φ2	0.03	Q-68	I-68	G-68
	45*12.5*22.5	20	φ2	0.06	Q-69	I-69	G-69
	45*12.5*12.5	10	φ1.5	0.018	Q-70	I-70	G-70
 Self masking continuous flowthrough cell	45*12.5*7.5	5	4*12	0.24	Q-73	I-73	G-73
	45*12.5*12.5	10	4*12	0.48	Q-74	I-74	G-74
	45*12.5*22.5	20	4*12	1	Q-75	I-75	G-75
	45*12.5*32.5	30	4*12	1.5	Q-76	I-76	G-76
	45*12.5*7.5	5	2*12	0.12	Q-79	I-79	G-79
	45*12.5*12.5	10	2*12	0.24	Q-80	I-80	G-80