





Cell-07							
Picture and Name	Outside Dim H*W*D(mm)	Path Length (mm)	Inside Width (mm)	Volume (ml)	Product Number		
					ES Quartz Glass 190nm-2500nm	IR Quartz Glass 250nm-3500nm	Optical Glass 340nm-2500nm
 Micro cell with frosted walls and with teflon stopper	45*12.5*7.5	5	2	0.35	Q-283	I-283	G-283
	45*12.5*12.5	10	2	0.7	Q-284	I-284	G-284
	45*12.5*22.5	20	2	1.4	Q-285	I-285	G-285
	45*12.5*32.5	30	2	2.1	Q-286	I-286	G-286
	45*12.5*42.5	40	2	2.8	Q-287	I-287	G-287
	45*12.5*52.5	50	2	3.5	Q-288	I-288	G-288
	45*12.5*102.5	100	2	7	Q-289	I-289	G-289
 Semi-Micro cell with frosted walls and with teflon stopper	45*12.5*7.5	5	4	0.7	Q-293	I-293	G-293
	45*12.5*12.5	10	4	1.4	Q-294	I-294	G-294
	45*12.5*22.5	20	4	2.8	Q-295	I-295	G-295
	45*12.5*32.5	30	4	4.2	Q-296	I-296	G-296
	45*12.5*42.5	40	4	5.6	Q-297	I-297	G-297
	45*12.5*52.5	50	4	7	Q-298	I-298	G-298
	45*12.5*102.5	100	4	14	Q-299	I-299	G-299
 Large cell	40*20*10.4	6.35	16	3.5	Q-321	I-321	G-321
	40*20*16.7	12.7	16	7.1	Q-322	I-322	G-322
	40*20*29.4	25.4	16	14.2	Q-323	I-323	G-323
	40*20*54.8	50.8	16	28.4	Q-324	I-324	G-324
	40*20*137.4	133.35	16	74.6	Q-325	I-325	G-325
	40*20*14	10	16	5.6	Q-326	I-326	G-326
	40*20*104	100	16	56	Q-327	I-327	G-327
 Micro cell with frosted walls and with teflon stopper	45*12.5*7.5	5	1	0.17	Q-373	I-373	G-373
	45*12.5*12.5	10	1	0.35	Q-374	I-374	G-374
	45*12.5*22.5	20	1	0.7	Q-375	I-375	G-375
	45*12.5*32.5	30	1	1	Q-376	I-376	G-376
	45*12.5*42.5	40	1	1.4	Q-377	I-377	G-377
	45*12.5*52.5	50	1	1.75	Q-378	I-378	G-378
	45*12.5*102.5	100	1	3.5	Q-379	I-379	G-379