

Cell-08							
Picture and Name	Outside Dim H*W*D(mm)	Path Length	Inside Width	Volume	Product Number		
		(mm)	(mm)	(ml)	ES Quartz Glass	IR Quartz Glass	Optical Glass
					190nm-2500nm	250nm-3500nm	340nm-2500nm
 Micro cell with frosted walls and with teflon stopper	45*12.5*7.5	5	3	0.5	Q-383	I-383	G-383
	45*12.5*12.5	10	3	1	Q-384	I-384	G-384
	45*12.5*22.5	20	3	2.1	Q-385	I-385	G-385
	45*12.5*32.5	30	3	3.15	Q-386	I-386	G-386
	45*12.5*42.5	40	3	4.2	Q-387	I-387	G-387
	45*12.5*52.5	50	3	5.25	Q-388	I-388	G-388
	45*12.5*102.5	100	3	10.5	Q-389	I-389	G-389
 Semi-Micro cell with frosted walls and with teflon stopper	45*12.5*7.5	5	5	0.85	Q-393	I-393	G-393
	45*12.5*12.5	10	5	1.7	Q-394	I-394	G-394
	45*12.5*22.5	20	5	3.5	Q-395	I-395	G-395
	45*12.5*32.5	30	5	5.2	Q-396	I-396	G-396
	45*12.5*42.5	40	5	7	Q-397	I-397	G-397
	45*12.5*52.5	50	5	8.75	Q-398	I-398	G-398
	45*12.5*102.5	100	5	17.5	Q-399	I-399	G-399
 Standard cell with lid and level bottom	45*12.5*3.5	1	10	0.35	Q-401	I-401	G-401
	45 x12.5*4.5	2	10	0.7	Q-402	I-402	G-402
	45*12.5*7.5	5	10	0.17	Q-403	I-403	G-403
	45*12.5*12.5	10	10	3.5	Q-404	I-404	G-404
	45*12.5*22.5	20	10	7	Q-405	I-405	G-405
	45*12.5*32.5	30	10	10.5	Q-406	I-406	G-406
	45*12.5*42.5	40	10	14	Q-407	I-407	G-407
	45*12.5*52.5	50	10	17.5	Q-408	I-408	G-408
	45*12.5*102.5	100	10	35	Q-409	I-409	G-409
 Micro cell with black walls and with lid	25*12.5*12.5	10	1	0.15	Q-424	I-424	G-424
	25*12.5*12.5	10	2	0.3	Q-434	I-434	G-434
	25*12.5*12.5	10	3	0.45	Q-444	I-444	G-444
	25*12.5*12.5	10	4	0.6	Q-454	I-454	G-454