





Cell-11							
Picture and Name	Outside Dim H*W*D(mm)	Path Length	Inside Width	Volume	Product Number		
		(mm)	(mm)	(ml)	ES Quartz Glass 190nm-2500nm	IR Quartz Glass 250nm-3500nm	Optical Glass 340nm-2500nm
 Micro cell with black walls	45*12.5*12.5	10	2	0.1	Q-674	I-674	G-674
 Absorptiometer rectangular	45*28*6.5	2.5	25	2.2	Q-682	I-682	G-682
	445*28*9	5	25	4.4	Q-683	I-683	G-683
	45*28*14	10	25	8.8	Q-684	I-684	G-684
	45*28*24	20	25	17.5	Q-685	I-685	G-685
	45*28*44	40	25	35	Q-686	I-686	G-686
	45*28*54	50	25	43.8	Q-687	I-687	G-687
	45*28*104	100	25	87.5	Q-688	I-688	G-688
	45*28*204	200	25	175	Q-689	I-686	G-689
 Rectangular with graded seal	45*12.5*12.5	10	10	3.5	Q-704	I-704	G-704
	Tube OD 8mm						
	Overall height 120mm						
 Fluorimeter rectangular with graded seal	45*12.5*12.5	10	10	3.5	Q-714	I-714	G-714
	Tube OD 8mm						
	Overall height 120mm						