


CM-1 CLARITY TESTER



CM-1 CLARITY TESTER

It tests clarify of injection and bottled medicine liquid.

Applicable standards

National standard: (Chinese Codex) /2015

Health department standard: WB-362 (B-121)-91

Corporate standard(Clarify tester): Q/12XQ0193-2010

Features:

The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor, and shading.

Flash eliminated. Illumination is adjustable.

Time and alarm automatically. Test continuously.

Used by one person.

Specifications:

Lamp three-primary-color fluorescent lamp

Illumination range : (1000~4000)LUX

Time range: (1~99)s

Power: 220V/ 50Hz/ 30W

Dimension: (660*180*500)mm³