

# CM-2 CLARITY TESTER



## CM-2 CLARITY TESTER

It tests clarify of injection and bottled medicine liquid.

Applicable standards

National standard: ( Chinese Codex ) /2015

Health department standard : WB-362 (B-121)-91

Corporate standard(Clarify tester ): Q/12XQ0193-2010

### Features:

The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor, and shading.

Flash eliminated. Illumination is adjustable.

Time and alarm automatically. Test continuously.

Used by two people at the same time.

### Specifications:

Lamp three-primary-color fluorescent lamp

Illumination range: (1000~4000)LUX

Time range: (1~99)s

Power: 220V/ 50Hz/ 50W

Dimension: (660\*380\*500)mm<sup>3</sup>