Cj-882A Magnetic Stirrer



Basic Info

Model NO.: CJ-882A

• Large Capacity: Digital Display

• Transport Package: 40L*38W*22h (Cm) G. W/N. W: 8kg/7kg

• Origin: China

• Material: Stainless Steel

• Trademark: JY

Specification: unit packingHS Code: 8479820090

Product Description

High Power Magnetic Stirrer:

1. Compact, large capacity, LED, safe & easy to use

2. CE, ISO13485/14001/9001

3. SGS Audited Factory Since 1984

Specifications of CJ-882A Thermostatic High Power Magnetic Stirrer	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	0-2600rpm
Stirring power	250W
Max. stirring capacity	10000ml
Workarea	30cm*30cm/20cm*20cm
Temp range	Room temp to 100 Celsius
Timer	120m
Digital display	yes
Large capacity, easy to use, cost efficient, widely applicable.	