

## Closed Heating Furnace/Multifunctional Furnace



### Characteristic

- Cold rolled steel shell
- Step-less temperature regulation
- Model DL: Fully closed heating tray, no open fire, and waterproof function.
- Fully closed slant control panel keeps its control system from corrosion of splashed and spilled liquid.

Technical parameters						
Model		DK-98-II				
Product specification	DL-I-15	single-joint(1kw)	single-joint(2kw)	double-joint	four-joint	six-joint
Heating type	heat transfer					
Shell material	spraying cold rolled steel					
Heating device	Teflon fully closed heating tray	resistance wire				
Rated power	1.5kw	1.0kw	2.0kw	1kw*2	1kw*4	1kw*6
Temp. control	rheostat step-less regulation					

Heating size(W*Dmm)	φ180	125*125	155*155	125*125	125*125	125*125
working size(W*D*H mm)	250*290*80	150*180*140	185*220*160	305*180*140	610*180*140	910*180*140
exterior size(W*D*H mm)	280*300*90	200*215*220	230*250*220	335*215*220	660*215*220	950*215*220
packing size(W*D*H mm)	AC220V/6.8A	AC220V/4.5A	AC220V/9.1A	AC220V/9.A	AC220V/18.2A	AC220V/27.3A
NW/GW	1.8/2kg	1.5/2kg	2.5/3.5kg	3/4kg	5.3/6.8kg	7.5/9kg