





Connecting Lens for Digital Camera and Video Camera



Connecting Lens for Digital Camera

Model		Connecting Way to Digital Camera	Connecting Size to Microscope Tube	Magnification	Suitable Digital Camera
NDPL-1		M42mm × 0.75	φ23.2mm φ30mm φ30.5mm	2×	Nikon SLR Camera
NDPL-2		M42mm × 0.75		2×	Canon SLR Camera
NDPL-3		M45mm × 0.75		7×	Nikon 5400
NDPL-4		M52mm × 0.75		7×	Canon A95
NDPL-5		M44mm × 0.75		7×	Ricoh GX8
NDPL-6		M52mm × 0.76		7×	Canon A620/ 630
NDPL-7		M41mm × 0.5		7×	Olympus 4040
NDPL-8		M28mm × 0.75		7×	Nikon 4500
NDPL-9		M52mm × 0.75		7×	Canon A650
NDPL-20		M48.5mm × 0.75		7×	Olympus C5060
NDPL-21		M52mm × 0.75		4×	Canon S5IS
NDPL-.....	Special Design for other models				
	T adaptor ring to connect the NDPL-1 and Nikon SLR camera				
	T adaptor ring to connect the NDPL-1 and Canon SLR camera				
	Connector for standard C mount adapter				
	φ30mm or φ30.5mm adaptor for different size of microscope tube				

Connecting Lens for Digital and Video Camera

Model		Connecting Way to Digital Camera	Connecting Size to Microscope Tube	Magnification	Suitable Video Camera
NDPL-10		C Mount	$\phi 23.2\text{mm}$ $\phi 30\text{mm}$ $\phi 30.5\text{mm}$	0.45×	DCE, MDCE, HDCE
NDPL-11		C Mount		0.35×	DCE-2, MDCE-5C HDCE-10C, HDCE-20C (Optional)
NDPL-12		C Mount		0.52×	HDCE-30C (Optional)
NDPL-13		C Mount		0.28×	DCE-2 (Optional)