

DB-IV Laboratory Heater/Laboratory Hot Plate/ Heating Plate

Basic Info

- Model NO.: DB-IV
- Type: Hot Plate
- Specification: ISO13485/14001/9001, SGS audited factory
- HS Code: 8516609000

- Application: Lab
- Trademark: JY
- Origin: Changzhou, China



Product Description

1. SUS thermostatic laboratory hot plate and heating mantle, uniform heating, energy saving, safe-to-use, non-deforming top plate, anti-corrosion. Various models available.
2. CE, ISO13485/14001/9001 Certified
3. As a manufacturer, we can offer most reasonable prices
4. More than over 20 year's history that we will provide products with high quality
5. Reliable and experienced delivery agent
6. We can develop according to customer's samples or drawings.

Continuous operating temperature of up to 300°C (572°F), with thermostat fitted, 400°C (752°F), with digital control fitted.

Heating provided via high-quality mineral insulated metal sheathed elements.
Fitted with a thermostat as standard (digital controller, optional extra)
Works off a single-phase power supply.

Top plate in mild steel (stainless steel or aluminum top plates are available as an optional extra)
Casing and terminal cover fabricated in stainless steel.
Larger or custom hotplates are available on request.

Specifications of DB-IV SUS Laboratory Hotplate	
Power Supply	AC 220V/110V 50Hz/60Hz
Digital Display	Yes
Working Area	L30 W40cm
Size	L30 W40 H12cm
Heating Power	1500W