

# DC-15B Super Thermostatic Water/Oil Bath



## Basic Info

- Model NO.: DC-15B
- Structure: Desktop
- Certification: ISO, SGS
- Trademark: JY/OEM
- Customized: Customized
- Material: Stainless Steel
- Application: School, Hospital, Lab
- Origin: China

## Product Description

DC water/oil baths are made of high-quality stainless steel, with characteristics of compact structure, microprocessor controlling, easy to operate, air-cooled compressor, overheat protection, overcurrent protection, and novel stainless steel circulator, making it widely applicable for oil, biochemical, electronic instrument, physics, engineering, medicine, foodstuff, and other applications.

Specifications:	
Power Supply	AC 220V/110V 50Hz/60Hz
Temperature control range	RT~200° C
Temperature fluctuation	± 0.5 Celsius
Working room size(mm³ )	280*250*200
Bath Depth(mm)	200
Circulator Flow(L/Min)	3
Overall Size(mm³ )	300*370*500
Digital display	Yes
Laboratory water bath of easy to use, quality stainless steel material, good reliability.	