DF-II Heat-gathering Magnetic Heating Stirrer



Basic Info

- Model NO.: DF-II
- Customized: Customized
- Structure: Desktop
- Type: Magnetic Stirrer
- Specification: SO13485/14001/9001, SGS Audited
- Application: Industry, School, Hospital, Lab
- Certification: ISO, SGS
- Material: Stainless Steel
- Trademark: JY
- Origin: Changzhou, China

Product Description

Magnetic Thermostatic Hotplate Stirrer/Electric Mixer

- 1. Compact, water/oil, safe, LED, thermostat
- 2. ISO13485/14001/9001, SGS Audited
- 3. Since 1984

Magnetic Hotplate Stirrer, an ideal essential device for laboratories, is widely applied in colleges and universities, environment protecting, researching, sanitation, anti-epidemic, petroleum, metallurgy, chemical and medical units.

Specifications of DF-2 Heat-gathering Magnetic Heating Stirrer	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	Startup to 2600rpm
Temperature range	Room temperature to 300 Celsius
Motor power	60W
Heating power	500W
Suitab	le for the water bath, oil bath.