

DHG-9031A Series



Oven (Natural Convection)

Features:

- * Stainless-steel chamber, semicircular arcs at corners for easy cleaning, and the space between the shelves in the chamber is adjustable.
- * With inspection window
- * LED display, PID controller with a timer.

Option:

- * Programmable controller: 15 periods 30 steps, 1 to 99hs 99mins for each period, fan speed 0 to 100% adjustable.
- * RS485 connector can connect computer and printer to record the parameters and the variations of temperature.
- * Independent temperature-limiting Alarm system.

Specification:

| Model | DHG-9031A | DHG-9051A | DHG-9091A |
|-------------------------|--------------------|------------------|------------------|
| Temperature Range | RT+10~200°C | | |
| Temperature Uniformity | 2.50% | | |
| Display Resolution | 0.1°C | | |
| Temperature Stability | ±0.5°C | | |
| Power Consumption | 600W | 900W | 1200W |
| Electrical Requirement | 220V 50HZ | | |
| Internal Dimension (mm) | 280*280*350 | 350*350*460 | 430*430*520 |
| Shelves | 27L | 56L | 96L |
| Shelves | 1(PCS) | 2(PCS) | 2(PCS) |
| Remark | natural convection | | |