

DJ-4D Electrical Stirrer



DJ-4D ELECTRICAL STIRRER

Features:

High automation. Easy operation.

Two muddlers (one small, one big) for different vessels and liquors.

Velocity data and time data can be preset at any time. the preset and real data can be display alternately in time.

Specifications:

Velocity range: (60~2000)rpm

Preset time range: (1~99)min

Maximal capacity: 15000ml (water)

Operational Continuity: >8HOURS

Power : 220V/ 50Hz/ 125W

Dimension: (360×240×720)mm³