

DZX Digital Display Thermo Vacuum Dry Box



Basic Info

- Model NO.: DZX
- Specification: CE, ISO13485/9001/14001 Certified
- HS Code: 8524300090
- Trademark: JY
- Origin: China

Product Description

Widely used for heating of substance of thermostat, easy labile and oxidizable. A certain amount of vacuum could be gained and filled with inert gas, ensuring quick drying for complex substances. Digital temp controller to set, display and control temperature.

Model:	6020	6021	6050	6051	6090	6210
Power:	~220V 50Hz					
Heating Power:	≥300W		≥1000W		≥1200W	≥2160W
ambient temperature :	5~40					
Temperature range :	50~180	50~150	50~250	50~200	50~250	
resolution:	0.1					

fluctuation:	±0.1					
Vacuum degrees:	≤133pa					
Matrial of working chamber:	1Gr18Ni9Ti	08F	1Gr18Ni9Ti	08F	1Gr18Ni9Ti	Gr18Ni9Ti
Size of working chamber:	300×300×275(mm)		415×345×370(mm)		450×450×450(mm)	560×640×600(mm)
Outer Size:	590×410×460(mm)		710×510×540(mm)		610×610×1400(mm)	720×800×1550(mm)
Load shelves:	2		2		2	3