

Db-III Laboratory Heater/Laboratory Hot Plate/ Heating Plate



Basic Info

- Model NO.: DB-III
- Trademark: JY
- Origin: China
- Stainless Steel Plate: Thermostat
- Specification: CE, ISO13485/14001/9001
- HS Code: 8419899090

Product Description

1. Bi-heat elements, temp sensor, uniform heating
2. As a manufacturer, we can offer most reasonable prices
3. More than over 20 year's history that we will provide products with high quality
4. Reliable and experienced delivery agent
5. We can develop according to customer's samples or drawings. Hot plates are generally used to heat glassware or its contents. Some hot plates also contain a magnetic stirrer, allowing the heated liquid to be stirred automatically.

Specifications of DB-III Thermostatic Hotplate for Laboratory	
Power Supply	AC 220V/110V 50Hz/60Hz
Temperature Control Range	Room Temperature to 300 Celsius
Digital Display	Yes
Working Area	L20 W30cm
Size	L20 W30 H12cm
Heating Power	1000W