

# Dhg-9101 Electrothermal Thermostatic Drying Oven, Air Dry Oven



## Basic Info

- Model NO.: DHG-9101
- Structure: Desktop
- Certification: ISO, SGS
- Type: Incubator
- Trademark: Gen Lab
- Origin: Changzhou, China
  
- Customized: Customized
- Material: Stainless Steel
- Application: School, Hospital, Lab
- Chamber: Stainless Steel
- Specification: 1pc in 1 box
- HS Code: 8419899090

## Product Description

### Electrothermal Thermostatic Drying Oven

1. Microprocessor and DC fan, varies sizes of work room, with thermal circle system.
2. ISO13485/14001/9001 Certified, SGS Audited Factory
3. Since 1984

Electrothermal thermostatic drying oven is made of stainless steel, with double-layer tempered glass. It applies over-heat alarm and protection, timer, dual digital display, microprocessor controller and permanent axial fan. Intelligent and high temp models have himer, power leakage protection, and over-heat alarm in accordance with JB/T5520-91.

<b>Model</b>	<b>Power Supply</b>	<b>Temperature Range</b>	<b>Temperature Fluctuation</b>	<b>Power</b>	<b>Work Room</b>
DHG-9101 Electrothermal Thermostatic Drying Oven	AC 220V/110V 50Hz/60Hz	50- 300celsius	± 1celsius	500W	350*350*350mm
	AC 220V/110V 50Hz/60Hz	50- 300celsius	± 1celsius	1000W	450*350*450mm
	AC 220V/110V 50Hz/60Hz	50- 300celsius	± 1celsius	1500W	550*450*550mm
	AC 220V/110V 50Hz/60Hz	50- 300celsius	± 1celsius	2000W	600*500*750mm