## GL-3250 Series magnetic agitator

- \*Cast of the yoke of traditional mind set of domestic magnetic agitators to promote latest fully molding upgraded products in novel appearance.
- \* Novel Integrated shell: The integral shell is for med with anti-flaming enhanced type plastic, resisting heat, acid, alkali and organic solvents;
- \*Unique heating design: The die-casting aluminum alloy heating plate exteriorly pained with Teflon material is unique in heat conducting effect, high in coldness, heat and corrosion resistance.
- \* Double effects in heat insulating: The double-layer plate at the bottom of heating plate may fully raise heat efficiency and stop heat from conducting to main frame.
- \*Integrally upgrading: The integrally formed shell with top convexity design may effectively prevent carelessly over flowed fluid in agitating operation from flowing into inside the main frame;
- \* Visualized temperature indication: The unique digital circuit design may accurately select and control heating plate temperature or directly

## Specification:

Model	GL-3250A	GL-3250B	GL-3250C
Speed	100~1200RPM	100~1200RPM	100~1200RPM
Max imagitation Capacity	10L	10L	5L
Temperature	RT~300C	RT~190C	No
Protection Temp	350C	220C	220C
Temp. control mode	Analogic	Digital	No
Precision of Temp. Control		+/-2C	
ExternalPT-100Sensor	No	Yes	No
Forward-backward rotation	No	Yes	No
Diameter of heating plate	150mm	150mm	
Maxim Heating function	350W	350W	
Power	220V/50Hz	220V/50Hz	220V/50Hz
Overall dimension	198*225*115mm	198*225*115mm	198*225*115mm
Weight	2.4KG	2.4KG	2.0KG



GL-3250A