GZ-A Digital Display High Power Laboratory Electric Stirrer/Overhead Stirrer



Basic Info

Model NO.: GZ-AStructure: DesktopCertification: ISO, SGS

• Trademark: JY

Specification: 1 pc in 1 boxHS Code: 8479820090

Customized: CustomizedMaterial: Stainless Steel

• Application: School, Hospital, Lab

• Transport Package: Carton

• Origin: China

Product Description

Brief Introduction:

1. High power, 20L of capacity, 2-step speed regulate, durable

2. CE, ISO13485/14001/9001

3. SGS Audited Factory Since 1984

Specifications of GZ-A Digital Display High Power Electric Stirrer	
Power Supply	AC 220V 50Hz/60Hz
Speed	60-500rpm(stepI), 240-2000rpm(step II)
Motor power	200W
Max. stirring capacity	20Liter
Maximum moment of stirring shaft:	1850gcm
Applicable medium viscosity	0-10000mpas
Large capacity, easy to use, cost-efficient, digital display, widely applicable.	