

# Ghp Water-Isolation Incubator; Water-Jacket Design Incubator



## Basic Info

- Model NO.: GHP
- Advantage: Easy Operation
- Specification: SO13485/14001/9001 Certified, SGS Audited Factory
- HS Code: 8419899090
  
- Type: Incubator
- Trademark: Jinyi
- Origin: Changzhou, China

## Product Description

1. Water-isolation type
2. ISO13485/14001/9001 Certified, SGS Audited Factory
3. Since 1984

Water-isolation type incubator is a sound instrument for thermostatic culture in sectors and departments of scientific research, medicine, biology, chemical, and agriculture.

Model	Power Supply	Temperature Range	Temperature Accuracy	Power	Work Room
GHP-50	220V 50HZ	RT+5-65celsius	± 1celsius	250W	35× 35× 40cm
GHP-80	220V 50HZ	RT+5-65celsius	± 1celsius	300W	40× 40× 50cm
GHP-160	220V 50HZ	RT+5-65celsius	± 1celsius	400W	50× 50× 65cm
GHP-270	220V 50HZ	RT+5-65celsius	± 1celsius	600W	60× 60× 75cm

## Features

1. It is water-jacketed with a glass door in the chamber for easy observation. When the glass door is opened, breeze circulation and heating system will automatically stop, which ensures no disadvantage of convection.
2. Microprocessor PID controller. The light alarm indicates over-temperature, low water-jacket level or high water-jacket level. Heating will be ceased when low water-jacket level occurs.
3. Independent alarm system for temperature-limiting ensures experiments run safely and no accident would happen (Option).
4. Stainless steel chamber. The water-jacketed heating method ensures a uniform temperature and also keeps the temperature constant for a long time after power off.
5. Recorder and printer are options which can print or record set parameters and trace variation curves of temperature (4~20mA standard signal). (Option)