


## HD-2 THICKNESS TESTER



### HD-2 THICKNESS TESTER

HD series of instruments are used for detecting a thickness of tablet and capsule.

Applicable standard

Corporate standard(Thickness tester) Q/12XQ0194-2010

#### **Features:**

Used for measuring a thickness of peak and wall of capsule & diameter of tablet and pill.

High-accuracy test with the most precise, high-quality homemade centigrade meter.

#### **Specifications:**

Test range: (0~30)mm

Test accuracy:  $\pm 0.01$ mm

Type of capsule: 00,0,1,2,3,4 type capsule

Range of test tablet: 0~30mm

Dimension: (200×80×50)mm<sup>3</sup>