


HD-4 THICKNESS TESTER



HD-4 THICKNESS TESTER

HD series of instruments are used for detecting a thickness of tablet and capsule.

Applicable standard

Corporate standard(Thickness tester) Q/12XQ0194-2010

Features:

Used for measuring a thickness of peak and wall of capsule & diameter of tablet and pill.
High-accuracy test with the most precise, high-quality homemade centigrade meter.

Specifications:

Test range: (0~30)mm

Test accuracy: ± 0.01 mm

Type of capsule: 00,0,1,2,3,4 type capsule

Range of test tablet : 0~30mm

Dimension: (160×60×50)mm³