

HH-W600 Three-Purpose Thermostatic Water Bath



Basic Info

- Model NO.: HH-W600
- with Cap: Stainless Steel Inner
- Specification: L60 W30 H18cm
- HS Code: 8419899090
- Type: Heating Type
- Trademark: JY
- Origin: China

Product Description

HH-W620 Digital Lab Water Bath

Digital Lab Water Bath

1. SUS inner, capacity of 30L, LED, high power
2. CE, ISO13485/14001/9001
3. Since1984

| Specifications of HH-W600 Thermostatic Laboratory Water bath/tank | |
|---|---------------------------------|
| Power Supply | AC 220V/110V 50Hz/60Hz |
| Temperature control range | Room temperature to 100 Celsius |
| Temperature coherence | 0.5 Celsius |
| Fluctuation of oil temperature | 0.5 Celsius |
| Heating power | 1500W |
| Working room size | L60 W30 H18cm |
| LED display | Yes |
| 30L capacity, easy to use, rich function configuration, high quality, wide applicability. OEM service is offered. | |