

HJ-4B Four-Unit Magnetic Laboratory Hotplate Stirrer



Basic Info

- Model NO.: HJ-4B
- Structure: Desktop
- Certification: ISO, SGS
- Type: Magnetic Stirrer
- Transport Package: Paper Packing
- Origin: Changzhou, China

- Customized: Customized
- Material: Steel
- Application: School, Hospital, Lab
- Trademark: JY
- Specification: 1 pc in 1 box

Product Description

Magnetic Thermostatic Hotplate Stirrer/Electric Mixer

1. 4 units, thermostatic, easy to use, 6 units are available
2. ISO13485/14001/9001, SGS Audited
3. Since 1984

Magnetic Hotplate Stirrer, an ideal essential device for laboratories, is widely applied in colleges and universities, environment protecting, researching, sanitation, antilipidemic, petroleum, metallurgy, chemical and medical units.

Specifications of HJ-4B Multi-unit Thermostatic Magnetic Hotplate Stirrer/Electric Mixer	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	0-2400rpm
Motor power	25W*4
Heating power	150W*4
Timer	120m
Temperature control range: room temp to 100° C	

