

HJ-5 Multifunctional Laboratory Magnetic and Electric Stirrer



Basic Info

- Model NO.: HJ-5
- Structure: Desktop
- Type: Magnetic Stirrer
- Specification: ISO13485/14001/9001, SGS Audited
- Customized: Customized
- Material: Steel
- Trademark: JY
- Origin: Changzhou, China

Product Description

Multifunctional Laboratory Magnetic and Electric Stirrer

1. Both magnetic and electric stirrer, easy to use
2. ISO13485/14001/9001, SGS Audited
3. Since 1984

Multifunctional Stirrer, an ideal essential device for laboratories, is widely applied in colleges and universities, environment protecting, researching, sanitation, antiepidemic, petroleum, metallurgy, chemical and medical units.

Specifications of HJ-5 Multifunctional Stirrer/Electric Mixer	
Power Supply	AC 220V/110V 50Hz/60Hz
Magnetic stirring speed	0-2400rpm
Electric stirring speed	0-2000rpm
Magnetic motor power	15W
Electric motor power	150W
Heating power	≤ 200W
Timer	No

