

HN-AT5 Low Speed Large Capacity Centrifuge



Bench top, large capacity, horizontal rotor, brush less motor, multi-capacity, sound structure, durable a large range of applicability’s.

HN-AT5 Low Speed Large Capacity Centrifuge

Technical data:

- Speed: start-4000 rpm (touch key control)
- Capacity: 100ml×4, 50ml×8, 15ml×24 (3 rotors)
- Max RCF: 2700 xg
- Power: 250W
- Timer: 0-99min
- Power Supply: 220V/110V 50Hz/60hz