

# HN-S12C High Speed Micro Hematocrit Centrifuge



HN-S12C

## Features:

1. As a manufacturer, we can offer most reasonable prices
2. More than over 20 year's history that we will provide products with high quality
3. Reliable and experienced delivery agent
4. We can develop according to customer's samples or drawings.

Fine brushless motor, easy& safe-to-use, low noise, touch key interface.

## Specifications of HN-S12C Benchtop Digital Micro hematocrit Centrifuge:

Power Supply	AC 220V/110V 50Hz/60Hz
Max speed	12000rpm Step less regulation
Max rcf	15300xg
Capacity	dia 1.5mm*75mm (24capillary tubes)
Timer	0-99min
Digital display	Yes, speed indication
Safety control	Automatic stop when lid is open during operation