## **HN-S802B**



## HN-S802B (80-2B) Low Speed Centrifuge

It is possessed of these special feature following: it is beautiful in appearance, large capacity but small in bulk and has a complete function and a stable performance, you can chose the speed that you need at the set range, and it can keep balance automatically, in addition, it has a little rate of increase in temperature and can be used efficiently and widely, It is a dependable instrument for qualitative analysis of serum, plasma and immunity in hospitals or chemical and biochemical laboratories.

## Specification:

1. Power Supply: AC 220V 50Hz (110V 60HZ)

2. Speed: 0-4000 rpm

3. Timer range: 0~60min

4. Max. RCF: 2325xg

5. Capacity: 12\*20ml

6. Motor power: 40 W

7. Weight: 9.5Kg