

HN-S803 (80-3) Centrifuge (LED)



Product Description

It has the advantages of small cubage, light weight, large capacity, low noise and so on. It can be operated conveniently. It is an ideal instrument for hospitals, labs to make qualitative analysis to serum, plasma, radio-immunity.

What we provide:

1. all the items to be of good quality.
2. all colors, designs, specification to be same as what the customer has confirmed with the fty
3. brand name under customers design will put on panel as contracted.
4. all parts to be intact and without any breakage on the goods.
5. the packing of all items to be of good quality, clean and attractive. make sure all the machine is well protected.
6. every machine will have a carefully inspection before packing and delivering
7. in consideration of there is no inspector from customer's side, we will take related photo and show customer the process of their order.

Specifications of HN-S803 Lab Digital Centrifuge

Power Supply: AC 220V/110V 50Hz/60Hz
Capacity: 20mlx6
Max rcf: 1430Xg
Timer: 0-60min
Max speed: 4000rpm
Digital display: Yes

User friendly, high quality motor drive, low noise, large capacity, steady performance, widely applicable for labs, biology, chemistry, physical, hospital, food and environmental protection fields.