

HY-4



HY-4 Multi-use Laboratory Shaker/Flask Shaker

1. SUS, high capacity, orbital
2. ISO13485/14001/9001 audited factory
3. Since 1984

The shaker is one with more advanced technique and reliable quality in its like. It has advantages of stable operation, low noise, stepless adjustable shaking speed and wide adjustable range. The spring multi-purpose jigger can match with diversified flasks to meet needs from different customers. It is an optimal equipment in research and production in biochemical industry, biological engineering, teaching and studying, microbiology and medical science and so on.

Easy to operate, with competitive price and high quality.

Specification:

Power Supply AC: 220V/110V 50Hz/60Hz

shaking amplitude: 20mm

shaking frequency : 0~300 times/min

Timer : 0-120min

Working load: 30kg

Working method : Reciprocating

Oscillation tray : 400*260mm