

# HY-PHS-25



## HY-PHS-25 PH METER

### Specifications :

Measuring Range pH: 0~14.00pH mV:0~±1400mV T:0~60°C

Accuracy pH: 0.05pH mV:±2%

Resolution: 0.05pH, 1mV

Manual Temperature Compensation Range 0~60°C

Input Impedance  $\geq 10^{12}\Omega$

### Using Condition :

Input Voltage : DC9V

Relative Humidity  $\leq 85\%$

Ambient Temperature : 5~40°C

No magnetic field interference, no significantly source

### Characteristics :

Automatic calibration

Automatic identification of standard buffer solution

Large screen LCD display