

HY-PHS-3C



HY-PHS-3C PH METER

Specifications :

Measuring Range pH: 0~14.00pH mV:0~±1999mV T:0~60°C

Accuracy pH: 0.01pH mV:±2%

Resolution: 0.01pH, 1mV

Manual Temperature Compensation Range: 0~60°C

Input Impedance $\geq 10^{12}\Omega$

Using Condition :

Input Voltage : DC 9V

Relative Humidity $\leq 85\%$

Ambient Temperature 5~40°C

No magnetic field interference, no significantly source

Characteristics :

Automatic calibration

Automatic identification of standard buffer solution

Large screen LCD display