

HZ-300



HZ-300 Rotary Shaker, featuring with compact structure and design, low-noise and smooth operation, is a laboratory instrument widely used in sections including biology, microorganism, biochemistry, virology, medical, environment protection and other scientific and research organizations.

Specification:

Power: AC 220V \pm 10%, 50Hz \pm 1Hz

Shaking Frequency: 40~300RPM/MIN

Amplitude: 20mm

Shaking Method: Rotary

Motor Power: 120W

Timer: 0~999min

Capacity: 250ml*20, 500ml*12, 1000ml*6

Size : 50cm \times 38cm \times 22cm

G.W/N.W: 26/33kg