

# Hh-S Digital Display Thermostatic Lab Oil Bath, Oil Bath



## Basic Info

- Model NO.: HH-S
- Trademark: JY
- Origin: Changzhou, China
- Type: Heating Type
- Specification: 1PC in 1 carton
- HS Code: 8419899090

## Product Description

Laboratory Oil Bath

1. S steel inner, LED, Thermostatic
2. ISO13485/14001/9001 Certified
3. Since 1984

| Specifications of HH-S Thermostatic Lab Digital Oil Bath   |                                 |
|--|---------------------------------|
| Power Supply   | AC 220V/110V 50Hz/60Hz          |
| Temperature control range  | Room temperature to 300 celsius |
| Temperature coherence  | 1Celsius                        |
| Fluctuation of oil temperature   | 2Celsius                        |
| Heating power  | 2000W                           |
| Working room   | 28X28X30(CM)                    |
| Digital display  | Yes                             |
| Easy to use, rich function configuration, high quality, wide applicability. Stainless steel made water tank. OEM service is offered. |                                 |