

# Hj-4A Thermostatic Magnetic Stirrer



## Basic Info

- Model NO.: HJ-4A
- Four Units: with Hotplate
- Specification: 1 pc in 1 box
- HS Code: 8479820090
- Type: Magnetic Stirrer
- Trademark: OEM/Unbrand
- Origin: China Mainland

## Product Description

Magnetic Hotplate Stirrer/Electric Mixer

1. 4 units, thermostatic, LED display speed and temp, easy to use, 6 units are available
2. ISO13485/14001/9001, SGS Audited
3. Since 1984Magnetic Hotplate Stirrer, an ideal essential device for laboratories is widely applied in colleges and universities, environment protecting, researching, sanitation, antilipidemic, petroleum, metallurgy, chemical and medical units.

Specifications of HJ-4A Multi-unit Thermostatic Magnetic Hotplate Stirrer/Electric Mixer	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	0-2400rpm, LED display
Motor power	25W*4
Heating power	150W*4
Timer	120m
Temperature control range: room temp to 100° C, LED display	